






MY Medikal

PRODUCT CATALOGUE

Update **Aug 2025**





Vision

We aim to become the leading-edge manufacturer, distributor and trusted partner for medical devices in the health care services and facilities.



Mission

1. To inspire healthier and well-being communities by providing easy sourcing healthcare products/medical devices.
2. To develop and deliver high-quality and safe healthcare products/medical devices worldwide by complying with the national and international standards and legal requirements.
3. To provide economically competitive safe healthcare products/medical devices to the market.
4. To keep customer satisfaction always at the forefront and at the highest level.



About us

My Medikal is a manufacturer and global distributor of various medical products such as syringes, gauzes, surgical and examination gloves both powdered and powder-free.

Our story begins in 2017 with the objective of presenting high-quality, safe and economically competitive products to the healthcare market.

Ideally located in Turkey, we have easy and strategic access to many countries in Asia and Europe which bring significant advantage when it comes to transportation and logistics.

Since the first day of its establishment, our company's value has brought together experienced and energetic people. As a result, we have constantly improved our processes and competencies with the young and dynamic management.

By the year of 2023, we are proudly operating with 2 factories in Izmir and Istanbul. Currently, My Medikal sells to nearly 40 countries around the world and with on-going plans of expansion.

mumu

Our Products

Sterile Surgical Gloves Powdered and Powder-free
(Latex & Non-Latex)

Latex Examination & Disposable Gloves
(Powdered & Powder-free)

Nitrile Examination Gloves Powder-free
Vinyl Gloves Powder-free

Syringe With & Without Needle (Luer Slip & Luer Lock)

Gauze Swabs & Lap Sponges (Sterile & Non-Sterile)

Non-Woven Gauze Swabs (Sterile & Non-Sterile)

mumu

Latex Powdered Disposable Gloves

FEATURES

Easy Donning
Outstanding Softness & Comfort
AQL 1.5
Smooth Surface
Made of Natural Rubber Latex



SPECIFICATION

Glove Details	Specification
Material	Natural Rubber Latex
Powder	> 2mg /glove
Surface	Smooth
Shape	Ambidextrous
Colour	Natural
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

PHYSICAL PROPERTIES

Glove Details	Properties
Weight (Size M)	4.8g (+/- 0.2g)
Freedom from Holes	AQL 1.5
Length (mm)	240 mm
Finger Thickness	Min. 0.08 (+/- 0.02 mm)
Palm Thickness	Min. 0.08 (+/- 0.02 mm)
Cuff Thickness	Min. 0.06 (+/- 0.02 mm)

Available Sizes & Palm Width (mm):

XS:70+/-10, **S:**80+/-10, **M:**95+/-10, **L:**111+/-10, **XL:**120+/-10

(According to ASTM D3578, ASTM D3767)

Packaging Style: 100 pieces according to weight or count.
20 inner boxes per master carton.

Country of Origin: MALAYSIA / THAILAND





mumu

Latex Powder-Free Disposable Glove



FEATURES

- Easy Donning
- Outstanding Softness & Comfort
- AQL 1.5
- Micro-Textured and Beaded Cuff
- Chlorinated

SPECIFICATION

PHYSICAL PROPERTIES

Glove Details	Properties
Weight (Size M)	4.8g (+- 0.2g)
Freedom from Holes	AQL 1.5
Length (mm)	240 mm
Finger Thickness (mm)	0.10 mm (+/- 0.02)
Palm Thickness (mm)	0.08 mm (+/- 0.02)
Cuff Thickness (mm)	0.06 mm (+/- 0.02)

Glove Details	Specification
Material	Natural Rubber Latex
Residual Powder	2mg /glove maximum
Surface	Textured
Shape	Ambidextrous
Colour	Natural
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

Available Sizes & Palm Width (mm):

XS:70+/-10, **S:**80+/-10, **M:**95+/-10, **L:**111+/-10, **XL:**120+/-10

(According to ASTM D3578, ASTM D3767)

Packaging Style: 100 pieces according to weight or count.
20 inner boxes per master carton



mumu

Latex Examination Glove Powderfree

FEATURES

High Level of Protection
Excellent Donning and Comfortable Fit
Chlorinated



SPECIFICATION

Glove Details	Specification
Material	Natural Rubber Latex
Residual Powder	Max. > 2mg /glove
Surface	Textured (Micro)
Shape	Ambidextrous
Colour	Natural / Pale Yellow
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

PHYSICAL PROPERTIES

Glove Details	Properties
Freedom from Holes	AQL 1.5
Weight (size M)	4.6g (+/-0.2g)
Length (mm)	240 mm
Finger Thickness (min)	0.08 mm
Palm Thickness (min)	0.08 mm
Cuff Thickness (min)	0.06 mm

SIZING

Size	Extra Small	Small	Medium	Large	Extra Large
Palm Width	76 +/- 3 mm	84 +/- 3 mm	94 +/- 3 mm	105 +/- 3 mm	113 +/- 3 mm

QUANTITIES

Packaging	Dimension	20 FT	40 HC	Packing
Outer carton	338 x 448 x 223 mm	850	2,004	20 Inner boxes / 2,000 pieces per carton
Inner Box	110 x 210 x 65 mm	17,000	40,080	100 pieces per inner box



MADE IN THAILAND



mumu Guard

Nitrile Examination and Protective Gloves

FEATURES

Latex-Free
Tactile Sensitivity
Chemical Resistance
Protection Against Viruses

QUALITY STANDARD:

Regulation (EU) 2017 / 7 45 PPE Regulation (EU) 2016/425 (EN ISO 374-1 : 2016 Cat 111)
EU Food Regulation 1935/2004 and 2023/2006 Regulation (EU) No. 10/2011 EN 455

PHYSICAL PROPERTIES

Glove Details	Properties
Elongation at Break - After Aging (%)	Min 400
Force at Break - Before & After Aging	≥ 6N
Freedom from Holes	AQL 1.5
Length (mm)	240 mm
Finger Thickness (mm)	0.09mm (+/- 0.02)
Palm Thickness (mm)	0.07mm (+/- 0.02)
Cuff Thickness (mm)	0.07mm (+/- 0.02)

SPECIFICATION

Glove Details	Specification
Material	Nitrile (NBR)
Surface	Fingertip Textured
Shape	Ambidextrous
Colour	Blue / Violet
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded



MADE IN CHINA

QUANTITIES

Packaging	Dimension	20 FT	40 HC	Packing
Outer carton	455 x 315 x 210 mm	960	2,267	20 Inner boxes / 2,000 pieces per carton
Inner Box	200 x 110 x 60 mm	19,200	45,340	100 pieces per inner box

mumu Guard

Nitrile Examination Glove Powderfree

FEATURES

Latex-Free
Tactile Sensitivity
Chemical Resistance
Protection Against Viruses



SPECIFICATION

Glove Details	Specification
Material	Nitrile (NBR)
Surface	Finger Textured
Shape	Ambidextrous
Colour	Black
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

PHYSICAL PROPERTIES

Glove Details	Properties
Force at Break - Before & After Aging	Min 6N
Elongation at Break - After Aging (%)	Min 400
Freedom from Holes	AQL 1.5
Length (mm)	Min 240 mm
Finger Thickness (mm)	0.09mm (+/-0.02)
Palm Thickness (mm)	0.07mm (+/-0.02)
Cuff Thickness (mm)	0.07mm (+/-0.02)

QUANTITIES

Packaging	Dimension	20 FT	40 HC	Packing
Outer carton	455 x 315 x 210 mm	960	2,267	20 Inner boxes / 2,000 pieces per carton
Inner Box	200 x 110 x 60 mm	19,200	45,340	100 pieces per inner box



MADE IN CHINA



mumu PROTECT

Nitrile Examination and Protective Powder-Free Gloves

FEATURES

PPE Regulation 2016/425, Cat III
MDR 2017/745, Class I

- This glove is suitable for food handling
- Black nitrile gloves are ideal for beauticians, mechanics, tattooist, factory workers and others.

PHYSICAL PROPERTIES

Glove Details	Properties
Donning Agent	Chlorinated
Elongation at Break - After Aging (%)	min. 400
Tensile Strength - Before & After Aging (MPA)	6 N
Freedom from Holes	AQL 1.5
Length (mm)	Min. 240
Cuff Thickness (mm)	0.09 +/- 0.03
Palm Thickness (mm)	0.10 +/- 0.03
Finger Thickness (mm)	0.13 +/- 0.03

SPECIFICATION

Glove Details	Specification
Material	Nitrile
Weight (Size M)	5.0g+- 0.3g
Surface	Finger Textured
Shape	Ambidextrous
Colour	Black
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

QUANTITIES / CONTAINER CAPACITY

Packaging Dimension	20 FT	40 FT	40 HC
315x250x245 mm	1150 Cartons	3000 Cartons	3600 Cartons
230x120x30mm	23,000 Boxes	60000 Boxes	72,000 Boxes

Packing Method:
50 pieces per box, 20 boxes per carton

Storage: Store in a cool dry place. Avoid excessive heat. Opened box should be shielded from exposure to direct sun or fluorescent lighting.
MADE IN CHINA

Remarks: Specification is subject to change without prior notice. Always confirm with our sales team.



NITRILE POWDER-FREE LONG CUFF BLUE EXAMINATION AND PROTECTIVE GLOVES

FEATURES

- MDR 2017/745 Class 1
- PPE Regulation 2016/425 Cat III
- European Commission Regulation (EU) No 2020/1245
- (EU) No 10/2011

SPECIFICATION

Glove Details	Specification
Material	Nitrile
Weight (Size M)	6.0 g \pm 0.3
Surface	Finger Textured
Shape	Ambidextrous
Colour	Blue
Sterilization	Non - Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

PHYSICAL PROPERTIES

Glove Details	Properties
Freedom from Holes	AQL 1.5
Length (mm)	Min. 280
Cuff Thickness (mm)	0.08 \pm 0.03
Palm Thickness (mm)	0.10 \pm 0.03
Finger Thickness (mm)	0.16 \pm 0.03

QUANTITIES / CONTAINER CAPACITY

Packaging	Dimension (mm)	Capacity 20FT	Capacity 40HC
Outer Carton	365 x 280 x 280 mm	1200 cartons	2400 cartons
Inner Box	270 x 135 x 35 mm	24000 boxes	48000 boxes

Packing Method: 50 pieces per inner box, 20 boxes per outer carton

Storage:

Store in a cool dry place. Avoid excessive heat. Opened box should be shielded from exposure to direct sun or fluorescent lighting.

MADE IN CHINA

mumu

Latex Surgical Gloves Powdered

FEATURES

MUMU Surgical Gloves designed with special features, elastic and very soft for long wear comfort. The micro-textured surface enhances slip-resistance giving better grip during long hour procedures.

SPECIFICATION

PHYSICAL PROPERTIES

Glove Details							Properties
Force at Break - Before & After Aging							≥ 9N
Freedom from Holes							AQL 0.65
Finger Thickness (mm)							0.17mm (+/-0.02)
Palm Thickness (mm)							0.15mm (+/-0.02)
Cuff Thickness (mm)							0.13mm (+/-0.02)
Size	6.0	6.5	7.0	7.5	8.0	8.5	
Length (mm)	260	260	270	270	270	280	
Width (mm)	77+/-5	83+/-5	89+/-5	95+/-5	102+/-5	108+/-51	

Weight. 17-18g/pair (Size 7.5)

QUANTITIES

Packaging	Dimension	20 GP	40 GP	40 HC
Outer carton	545 x 275 x 255 mm	800	1550	1800
Inner Box	131 x 131 x 240 mm	6.400	12.400	14.400
Inner Boxes	50 pairs per inner box			
Packing	8 Inner boxes / 400 pairs per carton			

Glove Details	Specification
Material	Natural Rubber Latex
Protein Content	32~48 ug/g
Surface	Micro - textured
Shape	Anatomic Shape (L/R), curved fingers
Colour	Natural
Sterilization	Irradiation
Shelf Life	5 Years
Cuff Finished	Beaded
Classification	
CE	Class IIa
Compliances	FDA 510K, CE mark, ISO 13485, ASTM D3577, EN 455:1:2:3.

Disclaimer: Product specification and artworks are subject to change without prior notice.

MADE IN MALAYSIA



mumu

Latex Surgical Gloves Powderfree

FEATURES

MUMU Surgical Gloves designed with special features, soft and elastic for long wear comfort. The micro-textured surface enhance slip-resistance for better grip during long hour procedures.

SPECIFICATION

Glove Details	Specification
Material	Natural Rubber Latex
Protein Content	32~48 ug/g
Surface	Micro - textured
Shape	Anatomic Shape (L/R), curved fingers
Colour	Natural
Sterilization	Irradiation
Shelf Life	5 Years
Cuff Finished	Beaded
Classification	
CE	Class IIa
Compliances	FDA 510K, CE mark, ISO 13485, ASTM D3577, EN 455:1:2:3.

QUANTITIES

Packaging	Dimension	20 GP	40 GP	40 HC
Outer carton	545 x 275 x 255 mm	800	1550	1800
Inner Box	131 x 131 x 240 mm	6.400	12.400	14.400
Inner Boxes	50 pairs per inner box			
Packing	8 Inner boxes / 400 pairs per carton			



PHYSICAL PROPERTIES

Glove Details		Properties					
Force at Break - Before & After Aging		≥ 9N					
Freedom from Holes		AQL 0.65					
Finger Thickness (mm)		0.19mm (+/-0.02)					
Palm Thickness (mm)		0.15mm (+/-0.02)					
Cuff Thickness (mm)		0.13mm (+/-0.02)					
Size	6.0	6.5	7.0	7.5	8.0	8.5	
Length (mm)	280	280	280	280	280	280	
Width (mm)	77+/-5	83+/-5	89+/-5	95+/-5	102+/-5	108+/-5	

Weight. 17-18g/pair (Size 7.5)

Disclaimer: Product specification and artworks are subject to change without prior notice.

MADE IN MALAYSIA





mumu

Latex Surgical Gloves Powdered

MUMU Surgical Gloves designed with special features, elastic and very soft for long wear comfort. The micro-textured surface enhances slip-resistance giving better grip during long hour procedures.

SPECIFICATION

Glove Details	Specification
Material	Natural Rubber Latex
Protein Content	40~50 ug/g
Donning Powder	Modified corn starch(USP)
Surface	Micro - textured
Shape	Anatomic Shape (L/R), curved fingers
Colour	Natural
Sterilization	Irradiation
Shelf Life	3 Years
Cuff Finished	Beaded

PHYSICAL PROPERTIES

Glove Details	Properties
Force at Break - Before & After Aging	≥ 9N
Freedom from Holes	AQL 0.65
Finger Thickness (mm)	0.20mm (+/-0.02)
Palm Thickness (mm)	0.18mm (+/-0.02)
Cuff Thickness (mm)	0.12mm (+/-0.02)

Size	6.0	6.5	7.0	7.5	8.0	8.5
Length (mm)	260	260	270	270	270	280
Width (mm)	77+/-5	83+/-5	89+/-5	95+/-5	102+/-5	108+/-51

Weight (average): 20g per pair (Size 7.5)

MADE IN CHINA

Inner Boxes	50 pairs per inner box
Packing	8 Inner boxes / 400 pairs per carton

Disclaimer: Product specification and artworks are subject to change without prior notice.

QUANTITIES

Packaging	Dimension	20 GP	40 GP	40 HC
Outer carton	535 x 265 x 410 mm	200,000	380,000	440,000
Inner Box	250 x 130 x 196 mm	4,000	7,600	8,800



mumu

Latex Surgical Gloves Powderfree

MUMU Surgical Gloves designed with special features, soft and elastic for long wear comfort. The micro-textured surface enhance slip-resistance for better grip during long hour procedures.

SPECIFICATION

Glove Details	Specification
Material	Natural Rubber Latex
Protein Content	32~48 ug/g
Surface	Micro - textured
Shape	Anatomic Shape (L/R), Curved fingers
Colour	Natural
Sterilization	Irradiation
Shelf Life	3 Years
Cuff Finished	Beaded
Classification	
CE	Class IIa
Compliances	FDA 510K, CE mark, ISO 13485, ASTM D3577, EN 455:1:2:3.

Disclaimer: Product specification and artworks are subject to change without prior notice.

QUANTITIES

Packaging	Dimension	20 GP	40 GP	40 HC	Inner Boxes	50 pairs per inner box
Outer carton	535 x 265 x 410 mm	200,000	380,000	440,000	Packing	8 Inner boxes / 400 pairs per carton
Inner Box	250 x 130 x 196 mm	4,000	7,600	8,800		



PHYSICAL PROPERTIES

Glove Details					Properties	
Force at Break - Before & After Aging					≥ 9N	
Freedom from Holes					AQL 0.65	
Finger Thickness (mm)					0.21mm (+/-0.02)	
Palm Thickness (mm)					0.20mm (+/-0.02)	
Cuff Thickness (mm)					0.16mm (+/-0.02)	
Size	6.0	6.5	7.0	7.5	8.0	8.5
Length (mm)	260	260	270	270	270	280
Width (mm)	77+/-5	83+/-5	89+/-5	95+/-5	102+/-5	108+/-5

Weight (average): 20g per pair (Size 7.5)

MADE IN CHINA





mumu

Synmax Vinyl Gloves

FEATURES

- Innovative formula of cost-efficient hand protection product.
- The performance of resistance to abrasion is ten times higher than the same nitrile gloves.
- Strong and durable performance
- No natural rubber latex protein to reduce risk of allergy

MANUFACTURING AND SAFETY STANDARDS

AQL	1
Quality Standards	Medical Grade
Powder Content	< 2mg/glove (according to EN 455-3)
Classification	Medical Device – Class I Regulation (EU) 2017/745 Personal Protective Equipment – Category III, Type B Regulation (EU) 2016/425
Compliances	EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009,
Viral Test	Test in accordance with ASTM F1671 & ISO 16604
Chemotherapy Drugs Permeation	Test in accordance with ASTM D6978
Chemical Substances Permeation	Test in accordance with EN 16523-1
Food Contact	Declaration of conformity for food contact in accordance with Regulation (EC) No 1935/2004 and with European Resolution Res AP (2004) 4, and Overall Migration Test in accordance with Commission Regulation (EU) No 10/2011
Shelf Life	5 years from the date of manufacturing

SPECIFICATION

PRODUCT DESCRIPTION	
Type of the glove	Single-use disposable gloves
Material	Vinyl
Not Made from Natural Rubber Latex	Yes
Shape	Ambidextrous
Colour	Blue
Powder Content	Powder-Free
Sterile	No
Cuff	Beaded

PHYSICAL PROPERTIES

Sizes	Weight	Length(mm)	Dimensions(mm)
S	4.0±0.3g	≥240	80 ±10
M	4.5±0.3g		95 ±10
L	5.0±0.3g		110 ±10
XL	5.5±0.3g		≥ 110
Thickness			
Finger	Palm	Cuff	Force at break
0.08±0.03mm	0.08±0.03mm	0.08±0.03mm	≥3.6N before and after aging

mumu

Clear Vinyl Gloves

Powdered & Powder Free



FEATURES

Latex free formulation, no natural rubber protein.
Chemical accelerator free.
Interior PU coated for easing donning.
Smooth finish for tactile sensitivity

SPECIFICATION

Glove Attributes	
Design	Ambidextrous, non sterile, beaded cuff
Donning Agent	PU coating on interior surface.
Colour	Clear, White, Yellow, Green, Red
Size	S, M, L, XL
Packing way	100pcs / box, 10 boxes / carton
Usage	Single use only

PHYSICAL DIMENSION

Width (+/- 5 mm)	S= 85 mm M = 95 mm L= 105 mm XL = 115 mm
Weight (+/- 0.3 gm)	Small = 3.9gm Medium = 4.0gm Large = 4.3gm X-Large = 4.5 gm

PHYSICAL PROPERTIES

Length (mm)	Minimum 230 mm
Thickness (mm)	
Cuff (mm)	N/A
Palm (mm)	Minimum 0.06 mm
Finger (mm)	Minimum 0.05 mm
Tensile Strength(MPa)	
Before aging	Minimum 9MPa
After aging	Minimum 9MPa
Elongation(%)	
Before aging	Minimum 300%
After aging	Minimum 300%



NATURAL RUBBER LATEX ORTHOPAEDIC

Surgical Gloves
Sterile Powderfree

PHYSICAL PROPERTIES

Glove Details						Properties
Force at Break - Before & After Aging						≥ 9N
Freedom from Holes						AQL 0.65
Finger Thickness (mm)						0.30mm (+/-0.03)
Palm Thickness (mm)						0.27mm (+/-0.03)
Cuff Thickness (mm)						0.20mm (+/-0.03)
Size	6.0	6.5	7.0	7.5	8.0	8.5
Length (mm)	280	280	285	285	290	290

PACKING INFORMATION

Outer carton	Long or square plastic peel pouch; Paper pouch		
Inner Box	Pouch	Inner box	Shipping Carton
	1 Pair / Pouch	50 Pair / Box	200 pairs or 400 Pairs / Carton

MADE IN CHINA

SPECIFICATION

Glove Details	Specification
Material	Natural Rubber Latex
Protein Content	32~48 ug/g
Surface	Micro textured, calibrated grip, medium
Shape	Anatomic Shape (L/R), curved fingers
Colour	Brown
Sterilization	Sterile, E-Beam / Gamma
Shelf Life	3 Years
Cuff Finished	Beaded cuff

Classification	
International Product Standards	EN 455 1-4, EN ISO 374 1-5:2016, EN ISO 21420:2019, ASTM D3577, ASTM D6978, ASTM F739, ISO 10993, AS/NZS 4179, ISO 10282, JIS T9107
Quality/ Environmental Standards	ISO 9001, ISO 13485:2016, ISO 14001

Medispo

Premium

POLYISOPRENE

SURGICAL GLOVES POWDERFREE

PRODUCT INFORMATION

Latex / Latex-free/ material	Latex-free, Polyisoprene
Color	White Cream
Cuff type	Beaded cuff
Texture / Grip	Micro textured, calibrated grip, medium
Sterile / Non-Sterile / Method	Sterile,E-Beam / Gamma
AQL value	0.65
International Product	EN 455 1-4, EN ISO 374 1-5:2016, EN ISO 21420:2019, ASTM D3577, ASTM D6978, ASTM F739, ISO 10993, AS/NZS 4179, ISO 10282, JIS T9107
Shelf Life	3 Years
Compliance to regulations	EU MDR 2017/745, USA 510k, MDSAP
PPE compliance / Type	PPE Regulation (EU) 2016/425, Type B

PACKING INFORMATION

Pouch material	Long or square plastic peel pouch; Paper pouch	
Packing Quantity	Pouch	Inner Box
	1Pair/Pouch	50 Pair/Box
	Shipping Carton	
	200 Pairs or 400 Pairs/Carton	

PRODUCT SPECIFICATION

Item No.	MPIG-5.5	MPIG-6.0	MPIG-6.5	MPIG-7.0	MPIG-7.5	MPIG-8.0	MPIG-8.5	MPIG-9.0
Size	5.5#	6.0#	6.5#	7.0#	7.5#	8.0#	8.5#	9.0#
Length (±5mm)	290	290	290	290	290	290	290	300
Thickness (±0.03mm)	Fingertip: 0.23		Palm: 0.23			Cuff: 0.16		

MADE IN CHINA



Syringe

With & Without Needles
STERILE, SINGLE USE, LATEX FREE



mumu

Syringe

With & Without Needles



Mumu Syringes are a sterile, single use device designed for the parenteral administration of pharmaceuticals. The intended purpose of device is to injection and aspiration of fluids from vials, vaccination applications, ampoules and parts of the body below the surface of the skin.

Class IIa Medicalpp Film+Paper Under Rule 2, Annex IX of Directive 93/42/EEC as amended by Directive 2007 / 47 /EC.

Needle Types		
18G 1"	18G x 1 1/2"	19G x 1 1/2"
20G x 1 1/2"	21G x 5/8"	21G x 2"
21G x 1 1/2"	22G x 1 1/4"	22G x 1 1/2"
23G x 1"	23G x 1 1/4"	25G x 5/8"
25G x 1"	26G x 1/2"	26G x 3/8"
27G x 1 1/2"	27G x 2"	30G x 1/2"

Syringe Details	Packaging	Shelf Box (Units)	Shipping Box (Units)
1 ml With & Without Needle	Medical pp film+paper	240	3600
2 ml With & Without Needle	Medical pp film+paper	350	2800
3 ml With & Without Needle	Medical pp film+paper	350	2800
5 ml With & Without Needle	Medical pp film+paper	250	2000
10 ml With & Without Needle	Medical pp film+paper	150	1200
20 ml With & Without Needle	Medical pp film+paper	100	800
50 ml With & Without Needle	Medical pp film+paper	40	320
60 ml With & Without Needle	Medical pp film+paper	40	320
50 ml Catheter	Medical pp film+paper	40	320
60 ml Catheter	Medical pp film+paper	40	320



mumu

GAUZE

STERILE & NON STERILE

Jumbo Gauze Rolls
Gauze Swabs
Lap Sponges



mumu

GAUZE

STERILE & NON STERILE



MUMU NON STERILE GAUZE SWABS

	Without Xray	With Xray
5 cm x 5 cm	8, 12 Ply	8, 12 Ply
7,5 cm x 7,5 cm	4, 8, 12, 16 Ply	8, 12, 16 Ply
10 cm x 10 cm	8, 12, 16 Ply	8, 12, 16 Ply
10 cm x 20 cm (+/- 2 cm)	8, 12 Ply	8, 12 Ply

MUMU NON STERILE LAPSPONGES

	Without Xray	With Xray
30 cm x 30 cm (+/- 2 cm)	8, 12, 16 Ply	8, 12, 16 Ply
40 cm x 40 cm (+/- 2 cm)	8, 12, 16 Ply	8, 12, 16 Ply
45 cm x 45 cm (+/- 2 cm)	8, 12, 16 Ply	8, 12, 16 Ply

MUMU NON STERILE GAUZE

Gauze With Xray	90 cm x 100 cm
Gauze Without Xray	90 cm x 100 cm

MUMU STERILE GAUZE SWABS - With Pouches

	Without Xray	With Xray
5 cm x 5 cm	8, 12 Ply	8, 12 Ply
7,5 cm x 7,5 cm	4, 8, 12, 16 Ply	8, 12, 16 Ply
10 cm x 10 cm	8, 12, 16 Ply	8, 12, 16 Ply
10 cm x 20 cm (+/- 2 cm)	8, 12 Ply	8, 12 Ply

MUMU NON STERILE LAPSPONGES - With Pouches

	Without Xray	With Xray
30 cm x 30 cm (+/- 2 cm)	8, 12, 16 Ply	8, 12, 16 Ply
40 cm x 40 cm (+/- 2 cm)	8, 12, 16 Ply	8, 12, 16 Ply
45 cm x 45 cm (+/- 2 cm)	8, 12, 16 Ply	8, 12, 16 Ply

Class IIa Medical PP Film+Paper Under Rule 2,
Annex IX of Directive 93/42/EEC as amended
by Directive 2007/47/EC.

Factory Management Certificates



Product Certificates/Standards



MDR 2017/745 Class I



MDD- 93/42/EEC -Class IIa





Disclaimer

Artworks or product specification are subject to change without earlier notice. Please consult and confirm with your sales representative.





MY Medikal

MY TİCARET VE MEDİKAL A.Ş.

Ömerli Mah. General Şükrü Koraltı Cad. No: 33
Arnavutköy / İstanbul / TÜRKİYE

T: +90 212 438 20 64 & F: +90 212 438 20 65
info@mymedikal.com.tr

www.mymedikal.com.tr